

# Simultaneous determination of eight drugs of abuse and codeine in saliva by liquid chromatography tandem mass spectrometry

K. Mortier, W. Lambert, J. Van Bocxlaer, D. Deforce,  
C. Van Peteghem and A. De Leenheer

Laboratory of Toxicology, Ghent University, Belgium



# Why saliva ?

## 1. Easy to sample

- Almost always available (? urine)
- No medical personnel needed (? blood)
- Painless, safer (? blood)
- Difficult to adulterate (? urine)

## 2. Saliva concentrations comparable to plasma ?

- Lipophilic basic drugs often ? C.
- Contamination ??

## 3. Relatively clean matrix

- Few lipids



# Measuring drugs of abuse in saliva

## 1. Often GC/MS

- Derivatization
- Extensive sample clean-up required

## 2. Often one group of substances

- Amphetamines
- Opiates
- Cocaine and benzoylecgonine

## 3. LC/MS seldom reported



# Why LC/ESI/MS QTOF?

## 1. LC/MS

- Simple sample prep
- No derivatization

## 2. Quadrupole Time of Flight

- Amphetamines
- Opiates
- Cocaine and benzoylecgonine

} In one run

## 3. Reliable identification

- Full spectrum without loss in sensitivity
- Accurate mass measurement



# Sample preparation (I)

## 1. Protein precipitation: easy and fast

- Suffers from matrix suppression (ESI)
  - ➔ better sample clean-up required
  - ➔ labelled internal standards ?
  - ➔ APCI (sensitivity?)

## 2. Mixed mode solid phase extraction

- Ideally suited for multiple drug extraction



# Sample preparation (II)

## Bond Elut Certify

200  $\mu$ l saliva + IS (+ standard) + phosphate buffer

Conditioning : 3 ml methanol  
2 ml phosphate buffer pH 6.0

Loading sample : 2 ml/min

Wash : 2 ml acetic acid (0.1M fast)  
: 2 ml methanol (fast)

Dry : 5 min

Elute : i-PrOH/CH<sub>2</sub>Cl<sub>2</sub>/NH<sub>4</sub>OH

Evaporate + HCl and redissolve in 200  $\mu$ l eluent



# Chromatography

Injection volume: 50  $\mu$ l

Column: phenyl column, 100 x 2.1 mm, 3  $\mu$ m

Flow: 0.2 ml/min, no splitting

Mobile phase:  $\text{NH}_4\text{Fo}$ , 10 mM, pH 5.0

Gradient: from 6 to 67.6% MeOH

Analysis time: 34 min.

Internal standards:

- Opiates : Butorphanol
- Amphetamines : MDMA
- Cocaine (and BE) : 2'-methylcocaine



# Mass spectrometry


| Compound         | Fragmentation                           | Cone V. | Coll. En |
|------------------|---|---------|----------|
| MDMPA            | 236 $\rightarrow$ 100 + 163             | 15      | 19       |
| Butorphanol      | 328 $\rightarrow$ 310                   | 37      | 32       |
| 2'-methylcocaine | 318 $\rightarrow$ 182 + 119             | 37      | 28       |
| Amphetamine      | 136 $\rightarrow$ 91 + 119              | 15      | 14       |
| Methamphetamine  | 150 $\rightarrow$ 91 + 119              | 15      | 18       |
| MDA              | 180 $\rightarrow$ 105 + 133 + 135 + 163 | 15      | 16       |
| MDMA             | 194 $\rightarrow$ 105 + 133 + 135 + 163 | 15      | 16       |
| MDEA             | 194 $\rightarrow$ 105 + 133 + 135 + 163 | 15      | 18       |
| Morphine         | 286 $\rightarrow$ 286                   | 38      | 13       |
| Codeine          | 300 $\rightarrow$ 300                   | 38      | 13       |
| Benzoyllecgonine | 290 $\rightarrow$ 168                   | 37      | 27       |
| Cocaine          | 304 $\rightarrow$ 182                   | 37      | 26       |






# Chromatogram I: Amphetamines

Blank saliva sample spiked at the 10 ng/ml level



# Chromatogram II: Opiates and cocaine

Blank saliva sample spiked at the 10 ng/ml level



# Validation (at low and high C)

Extraction recovery: Morphine 58 – 64%

Others: > 74%

No Matrix suppression: (5 diff. individuals)

Precision: Within and between day: < 16.8%

Accuracy: < 12%


LOQ: 2 ng/ml saliva

Calibration: Quadratic curves 2 – 100 ng/ml



# Samples

- 41 saliva samples from suspected drug users
- Samples: peak area of IS **irreproducible**
- Calibrators: peak area of IS **reproducible**



# Contamination origin

- Contamination from polymer releasing combinations of monomers
- Origin: Sample device (HSW Henke – Sass, WOLF GMBH)
- Suppression of later eluting compounds (IS)

## Conclusion

- Method validated for the determination of drugs of abuse in saliva, obtained by spitting
- Sampling device should be taken into account when developing analytical procedures

# Acknowledgements

- This work was supported by the EU by grant SMT4-CT-98-2257